

Patient Information from



URINE COLLECTION (24-HOUR) P-SPLITTER® DEVICE


This collection device is used when your physician has ordered more than one 24-hour urine test requiring different preservative conditions.

To avoid backpressure in bottles and ensure accurate splitting loosen bottle caps slightly to allow air pressure to be released when urine is poured into bottles.




Please **report to Meadville Medical Center Registration** before delivering the specimen to the lab. This is necessary for complete specimen identification and processing.

PROPER POURING INSTRUCTIONS:

RIGHT:

To assure equal splitting, be sure to "follow through" with a quick and uninterrupted pouring motion, to keep the urine level as high in the tube as possible.

WRONG:

If sample is poured slowly or "dribbled" into device, unequal splitting will occur.

NOTE: After sampling is complete, P-Splitter funnel should be removed from containers. Red plugs should then be inserted into container holes and covered thoroughly with attached labels. P-Splitter funnel should then be inserted into the enclosed plastic bag for disposal.